
# iVascular® iVolution

EN	INSTRUCTIONS FOR USE SELF-EXPANDING PERIPHERAL STENT SYSTEM	02
ES	INSTRUCCIONES DE USO SISTEMA DE STENT PERIFÉRICO AUTOEXPANDIBLE	03
PT	INSTRUÇÕES DE USO SISTEMA DE STENT PERIFÉRICO AUTO-EXPANSÍVEL	05
IT	ISTRUZIONI PER L'USO SISTEMA DI STENT PERIFERICO AUTOESPANDIBILE	07
FR	INSTRUCTIONS D'UTILISATION SYSTÈME D'ENDOPROTHÈSE PÉRIPHÉRIQUE AUTOEXPANSIBLE	09
DE	GEBRAUCHSHINWEISE PERIPHERES SELBST-EXPANDIERENDES STENTSYSTEM	11
RU	ИНСТРУКЦИИ ПО ИСПОЛЬЗОВАНИЮ - СИСТЕМА САМОРАСШИРЯЮЩЕГОСЯ ПЕРИФЕРИЧЕСКОГО СТЕНТА	13
CS	NÁVOD K POUŽITÍ SYSTÉM SAMOEXPANDIBILNÍHO PERIFERNÍHO STENTU	15
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ ΣΥΣΤΗΜΑ ΑΥΤΟΔΙΑΤΕΙΝΟΜΕΝΟΥ ΠΕΡΙΦΕΡΙΚΟΥ ΣΤΕΝΤ	16
PL	INSTRUKCJA OBSŁUG KOMPLEKT DLA USTANOWKI	18
NL	GEBRUIKSAANWIJZING ZELF-UITZETTEND PERIFERE STENTSYSTEEM	20
SE	ANVÄNDNINGSSINSTRUKTIONER SJÄLVEXPANDERANDE PERIFERT STENTSYSTEM	22
TU	KULLANMA TALİMATLARI KENDİNDEN GENİŞLEYEN PERİFERİK STENT SİSTEMİ	24

## ENGLISH

### 1. Product description

The self-expanding peripheral stent system from LVD Biotech SL is a self-expanding stent made of a nickel-titanium alloy, premounted on a coaxial delivery catheter.

The stent has a series of radiopaque markers at each end to mark its length and identify the position and location within the artery.

The stent is a flexible prosthesis that reaches its maximum diameter when removed from the delivery system, inside the artery. Once inside the artery, it adapts to it by exerting a constant radial force that maintains the vessel's permeability. The stent is made from a tube that is cut using a laser technique and then expanded to the required final diameter. The surface is then polished to achieve a smooth, shiny finish.

The stent's delivery system is a coaxial catheter compatible with a 0,035 guide. At the proximal end there is a handle, which when pushed, removes the stent's outer covering sheath releasing the stent.

The handle has a slow-delivery system, to begin delivery of the stent in an accurate manner, and a quick-delivery system which, once the stent is placed against the vessel's wall, can be removed quickly. The handle also contains a safety system to prevent the stent from being released or the sheath from moving before it is removed.

The system ends in a soft, atraumatic tip to avoid damaging the artery during its advance.

In the distal portion of the release system, there is a series of radiopaque markers which assist in the placement of the stent on the lesion and in the deployment stage of the same. These markers are:

- Two fixed markers which define the position of the stent, one situated distal to the stent and the other situated proximal to the stent
- A mobile marker, situated at the extreme distal end of the retractable sheath. Initially, this marker is aligned with the distal fixed marker, in such a way that, under the fluoroscope, only two markers in the system are observed: the distal and the proximal defining the stent. During the release of the stent, the marker moves backwards with the sheath: at this time, three markers are observed: the distal marker, the mobile marker that indicates the length of the released stent, and the proximal fixed marker. When the stent is completely released, the mobile marker of the retractable sheath is aligned with the proximal marker and thus, only two markers are seen.

The system comes in two lengths: 80 cm and 140 cm.

The self-expanding stent is available in different sizes, as summarised in the following table.

The variant selected for the procedure will depend on the diameter of the vessel to be treated. Depending on the variant, the specifications of the introducer and the guiding catheter to be used are defined.

### Contents

- Peripheral self-expanding stent mounted on its delivery system. Packed in a blister and bag to protect sterility.
- Implant card.
- A leaflet with instructions for use.

### 2. Indications

The peripheral stent system is indicated for the treatment of de novo or restenotic atherosclerotic lesions in peripheral arteries located under the aortic arch, that is, and for palliation of biliary tract malignant stenosis with a nominal diameter ranging from 4.5 to 9.5 mm.

### 3. Contraindications

- Patients with a contraindication for anticoagulant/antiplatelet therapy.
- Patients clearly diagnosed as allergic to heavy metals. (Specifically to nickel).
- Patients with hemorrhagic disorder.
- Presence of recent undissolved, disorganized thrombus or of embolic material.
- Renal failure or allergy to contrast media.
- Severe ascites.
- Threat of vital lateral branch occlusion.

### 4. Warnings

- Stents must only be implanted in hospitals equipped with emergency facilities and staff able to switch the operation to open surgery, or failing

this in hospitals with such services close at hand.

- The product must only be used by physicians well trained in angioplasty techniques, acquainted with percutaneous intervention and stent implantation.
- This device is designed and intended for single use in only one patient. DO NOT RSTERILIZE OR REUSE. Do not reprocess. Reutilization of the product in another patient may cause crossed contamination, infections or transmission of infectious diseases. Reutilization of the product in another patient may cause alterations thereof and limit its effectiveness.
- The product is supplied sterile. Check the expiry date and do not use products which have passed this date.
- Administer adequate medical therapy to the patient: anticoagulants, vasodilators, etc., in accordance with the intravascular catheter insertion protocol.
- Use aseptic techniques when removing the product from its primary packaging.
- Do not dry with gauze.
- Do not expose the delivery system to organic solvents.
- Select the right size in diameter and length according to the size of the lesion observed under fluoroscopy.
- Advance the product along the guidewire under fluoroscopy, until reach the treatment area. Do not allow the product to advance without the guidewire inside it.
- Do not superimpose stents of different compositions.
- Do not expose the product to hospital waste material.

### 5. Precautions

#### General Precautions

- Store in a dry, fresh place, away from direct sunlight.
- Inspect the pack before opening it. If it is found to be damaged, do not use the product.
- The minimum size of the introducer with which the system may be used is indicated on the product's label. Do not use introducers that are smaller than the size indicated. In addition, in case of long and/or braided introducer sheaths, a larger French-size than indicated on the label may be necessary to ensure success of the procedure.
- After use, this product may be a biological hazard. Handle and dispose of it in accordance with accepted medical practice and relevant local, state or federal laws and regulations

#### Precautions when preparing the product

- Inspect the product before use to make sure it is in good condition, and discard any device which does not inspire confidence. If you notice that the stent has begun to deploy or that the sheaths do not cover the stent, discard the device.
- The stent's length must coincide with the length of the lesion. Preferably, it should be a little longer so as to ensure that the lesion is totally covered.
- Do not attempt to remove the stent from its delivery system.
- Do not use the system if the safety device has moved or the deployment button has moved from its initial point.
- Do not unlock the handle before the stent has been placed in the precise location, because if it is unlocked before, the stent may be delivered in the wrong position.

#### Precautions during advance-positioning-deployment-withdrawal of the product

- Always advance the device under fluoroscopic visualization. If any resistance is met while advancing, stop and determine the cause before continuing. If the stent can not reach or pass through the lesion, withdraw the whole system as a single unit, including the introducer and the guidewire.
- If any resistance is met while advancing, stop and determine the cause before continuing.
- Before stent deployment, apply tension to the delivery system in order to avoid that the stent be implanted in an incorrect place.
- Once expansion of the stent has begun, the system cannot be withdrawn or moved.
- Do not handle, advance or withdraw neither the catheter nor the

guidewire while the stent has started its deployment.

- The stent is not designed to be reintroduced in the sheath once deployment has started.
- Once the stent has been placed against the vessel's wall, do not try to reposition or move the stent as it may drag and damage the artery.
- If resistance of any kind is noted while accessing the lesion, while extracting the deployment system or while removing a stent which could not be implanted, the whole system should be withdrawn as a single unit.
  - Do not retract the delivery system inside the introducer or the guiding catheter since the stent may be released into the vasculature.
  - Advance the guidewire forward in the anatomy, as far as possible.
  - Withdraw and extract the introducer (or the guiding catheter) and the stent system together, as a single unit
  - Finally, withdraw the guidewire, or in the event that you wish to begin implantation again keep the guidewire in place and place further devices.
- Positioning the stent in an important bifurcation may obstruct or hamper future access to lateral branches.
- For treatment of multiple lesions, start by the most distal lesions and continue onto the proximal ones.
- Immediately after an implant, if it is necessary to re-cross the stent, great care must be taken not to push it with the guides or balloons so as not to dislodge it.
- Corrosion by friction of superimposed prosthesis has not been demonstrated.

#### Precautions: MRI Safety[1]:

Thorough non-clinical testing security was proved by magnetic resonance (MR) imaging of nickel titanium alloy stents under certain conditions (MR conditional). The conditions of the MR test used to assess these devices were the following: for magnetic field interactions, a static magnetic field strength of 3.0 Tesla, with a maximum spatial gradient of 3,3 T/m; for MR-related heating, a maximum whole body averaged specific absorption rate (SAR) of 2,0W/kg for 15 minutes.

No clinical or non-clinical testing has been conducted to rule out the possibility of stent migration at field strengths of over 3.0 Tesla. MR imaging quality may be affected if the area of interest is adjacent or very close to the stent place.

#### 6. Possible adverse effects/complications

Possible adverse effects and/or complications which could arise before, during or after the procedure include the following:

- Death
- Infarction/angina
- Stroke / embolism / thrombosis
- Emergency bypass surgery
- Endocarditis
- Total occlusion of the artery / thrombosis / stent embolisation
- Perforation or dissection of the segment treated
- Air embolism
- Arteriovenous fistula
- Renal failure
- Migration of endoprosthesis / bad positioning / partial deployment
- Allergic reaction to the metals or contrast media
- Reocclusion of the segment treated: restenosis
- Arterial spasm
- Haemodynamic deterioration
- Arrhythmias
- Aneurysm or pseudoaneurysm
- Infections
- Hypotension/hypertension

[1] This data has been obtained bibliographically, searching for currently approved stents of equal composition. It is recommended that MRI is not performed unless strictly necessary for the patient. If the patient can wait, it is better to perform the MRI when the stent has endothelialized, around four weeks after implantation.

- Local haemorrhage with haematoma at the insertion site

Specifically for biliary use

- Abscess
- Additional intervention due to, but not limited to:
  - Endoprosthesis migration
  - Bad positioning
  - Partial deployment
- Biliary tract occlusion
- Perforation of biliary tract
- Cholangitis
- Pancreatitis
- Peritonitis
- Sepsis

## 7. How to use

### 7.1 Materials required

- Standard heparinised saline solution.
- Contrast medium (use contrast media indicated for intravascular use).
- Medication appropriate to the procedure.
- 0.035" Guidewire, do not use other guidewire sizes.
- Introducer with haemostasis valve in the size indicated on the label. If a long or braided introducer sheath is used, it is recommended that a larger French-size than the one indicated on the label be used.
- Several standard 10-20cc syringes of saline solution for flushing the system.
- Depending on user preference, it is possible to use a guiding catheter over the introducer, in which case it will be necessary to place a hemostatic valve in the guiding catheter's luer connector, and then continue inserting the device inside the catheter.

### 7.2 Preparing the System

- The diameter size of the stent has to fit to the diameter of the reference vessel (see the product specifications table in paragraph 1. Description, to decide the implanted stent diameter depending on the diameter of the vessel to be treated)
- The size of the stent should coincide in length with the length of the lesion. It is preferably slightly longer to ensure complete coverage of the lesion.
- Remove the catheter from the protection blister. Check that the stent is covered by the sheath. Do not remove the handle's safety system until deployment is going to begin.
- Flush or rinse the guide lumen with sterile saline injecting it with a syringe from the proximal connector to the handle, until it comes out the tip. Do not dry with gauze. Avoid touching the stent area.

DANGER: if the solution does not exit through the tip of the delivery catheter, the catheter is not in good condition and must not be used.

## 7.3 Implantation Procedure

### Preliminary steps

- Place the introducer in accordance with the manufacturer's instructions. If you wish to use the guiding catheter, insert it through the introducer following the manufacturer's instructions and place a haemostasis valve on its connector.
- Flush the guidewire with sterile heparinized saline. Do not dry with gauze.
- Introduce the guidewire through the introducer's haemostasis valve (or the guiding catheter), according to the manufacturer's instructions. Advance the guidewire throughout the vessel toward and through the lesion under fluoroscopic visualization.
- Before implanting a self-expanding stent, it is recommended that a balloon predilation is performed using PTA techniques to prepare the lesion and estimate its condition as well as the adequacy of the stent to be implanted. For balloon dilation, follow the manufacturer's instructions for use attached to the device.

### Implantation technique

- Insert the proximal end of the guidewire in the distal tip of the delivery

system until the guide exits through the proximal connector of the delivery system.

- Carefully insert the system through introducer under fluoroscopic visualization (if you are using a guiding catheter, keep the haemostasis valve open when the system passes through it). Once the introducer is through, make sure by fluoroscopic visualization that the sheath that covers the stent has not moved and the stent is still compressed inside it (movement of the sheath can be detected, if you see three radiopaque markers instead of two in the release system under fluoroscope, one of them could be observed on the stent, at a different height depending on the movement of the sheath, and indicating that the stent has begun to be released). If you notice that the stent is beginning to deploy, the system must be removed along with the introducer as a single unit according to the indications of the precautions section in this leaflet.
  - Carefully advance the stent system on the guidewire, with the help of fluoroscopy, until the area to be treated is reached. Make sure that the guidewire goes out through the tip of the delivery system.
  - Advance the stent beyond the lesion and pull back to tension the system and eliminate slacks. Preferably, the fixed sheath should be inserted into the introducer's haemostasis valve, if this is not the case, tighten the system with special caution. Make sure the stent is on the lesion to be treated by checking that the lesion is between the distal radiopaque marker of the tip aligned with the sheath marker and the radiopaque marker proximal to the stent, contained in the delivery system. Before expanding the stent, check by means of high resolution fluoroscopy that it has not been damaged or moved during advancement, otherwise, proceed to withdraw the system according to the withdrawal procedure indicated in the precaution section in this leaflet
  - Remove the security system to unlock the handle release system.
  - To release the stent can be used a slow speed and fast speed. Slow speed is driven by the wheel while the rapid deployment occurs when the button is pressed down and pulled back. We recommend starting the release of the stent using the thumbwheel on low speed until this is attached to the vessel wall. Once successful, you can use the fastest speed to expedite the procedure.
  - Hold on with one hand the proximal end of the handle with the thumb placed on the wheel. The other hand should hold the fixed sheath delivery system. Check that the introducer has not moved.
  - Start the removal of the cover sheath stent turning the dial to the right following the arrow. Before discovering the stent, the sheath will move about 3 mm, from there, the stent begins to release.
  - Before the stent is placed against the vessel's wall, it is possible to slightly reposition the system. Once it has been placed, do not move forward or backward, since the system may be damaged, or injury caused to the patient.
  - When you begin to observe that the marker of the sheath has withdrawn backward, the release of the stent has begun. The moving marker easily indicates the length of the released stent.
  - Once the stent has been placed against the vessel, and depending on the user's preferences, it is possible to begin the quick-delivery movement pressing down and moving back the release button. Avoid abrupt movements in this action since it may cause the stent to lengthen when the sheath is dragged on it.
  - Once it is confirmed that the stent is totally deployed, which as a guide, you can tell it is finished when the marker of the sheath is aligned with the proximal marker of the release system, remove the delivery system through the introducer
  - Perform an angiography (or cholangiography in biliary case) to confirm correct positioning of the stent.
  - Maintain the guidewire across the dilated stenosis for 30 minutes after the angioplasty. In case a thrombus is visualized around the stent, administer thrombolytic agents. When the angiography confirms dilation, carefully withdraw the guidewire.
  - Leave the introducer in situ until the haemodynamic profile is back to normal. Suture if necessary, in the usual manner.
- Further dilation of stented segments**
- If the stent size with respect the vessel is insufficient, the technique of stent post-dilation may be applied with a post-dilation balloon catheter. To this end, a balloon catheter is inserted in the stent area (it might be non-compliant) and the dilatation is performed according to the manufacturer's recommendation. Do not dilate the stent beyond its

maximum recommended artery implantation diameter.

### Withdrawing an unexpanded stent

- If the stent does not go across the lesion, it may be necessary to withdraw it unexpanded. To this end, follow the instructions recommended in the precautions in this leaflet and withdraw the system as a single unit.

## 8. Warranty

The device and all its components have been designed, manufactured, tested and packed with a maximum level of precaution. LVD Biotech warrants the device until its expiration date, provided the packaging is not broken, tampered with or damaged.

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© iVascular - cod. MP45012 - rev2 Nov 2016